



CRITIKON SENSA-CUF

Blood Pressure Cuffs

The CRITIKON™ SENSA-CUF™ is a highly versatile solution for a wide range of extended-care settings. Made using strong, breathable, easy-to-clean material, the SENSA-CUF is designed to maximize patient comfort in low-acuity areas in a hospital where risk of contamination are diminished, and durable enough for multi-patient use.



A trusted choice for clinical excellence

- Cuffs are held to the same quality standard as all GE patient monitors, with repeated inflations to ensure accurate readings
- Large bold print helps you quickly choose the right cuff, enabling a more accurate measurement

Constructed for patient safety and durability

- Soft, breathable, easy-to-clean material is not made with PVC, making them safer for the environment, or natural rubber latex, making them safer for your patients
- DINACLICK™ connector meets international standards and prevents misconnections to luer fittings

Designed to improve workflow

- Easy-to-use, fast DINACLICK connector enables nurses to link the cuff to the monitor faster than other dual-tube connectors
- Cuffs compliant with AAMI and AHA NIBP cuff sizing standards for prompt, correct sizing
- A full range of adapters and connectors are compatible with automated monitors and a variety of manual devices, including wall manometers, enabling a universal connection system throughout the facility

Assortment packs available

- For low-volume cuff utilizers, several types of assortment packs are available for easier ordering and more optimal inventory management



DINACLICK
Connector



SENSA-CUF			GE Healthcare			Non-GE Healthcare
						
Size	Color	Limb Circumference	2-Tube Screw Connector	2-Tube DINACLICK Connector	2-Tube Submin Connector	1-Tube Bayonet Connector
Infant	Orange	8-13 cm	2458	SEN-P1-2A	2482	2470
Child	Green	12-19 cm	2460	SEN-P2-2A	2484	2472
Small Adult	Lt. Blue	17-25 cm	2462	SEN-A1-2A	2486	2474
Small Adult Long	Lt. Blue	17-25 cm	2463	SEN-A1-2A-L	2487	2475
Adult	Navy	23-33 cm	2464	SEN-A2-2A	2488	2476
Adult Long	Navy	23-33 cm	2465	SEN-A2-2A-L	2489	2477
Large Adult	Rose	31-40 cm	2466	SEN-A3-2A	2490	2478
Large Adult Long	Rose	31-40 cm	2467	SEN-A3-2A-L	2491	2479
Thigh	Brown	38-50 cm	2468	SEN-T1-2A	2492	2480

Note: 5 cuffs/box

SENSA-CUF Assortment Pack	Part Number	Description
2-Tube Screw Connectors	2494	1 Small Adult, 1 Adult, 1 Large Adult (3 cuffs/box)
2-Tube DINACLICK Connectors	2059300-001	1 Small Adult, 1 Adult, 1 Large Adult (3 cuffs/box)

To learn more about GE Healthcare's comprehensive cuff solutions, visit us at cuffs.gehealthcare.com



Imagination at work

www.gehealthcare.com. Product may not be available in all countries and regions. Contact a GE Healthcare Representative for more information. Data subject to change.

© 2016 General Electric Company.

GE, GE Monogram, Imagination at work, CRITIKON, DINACLICK and SENSA-CUF are trademarks of General Electric Company.

Reproduction in any form is forbidden without prior written permission from GE. Nothing in this material should be used to diagnose or treat any disease or condition. Readers must consult a healthcare professional.